

Product Datasheet

Anti-Human IL1RAP Monoclonal Antibody, Azide Free Clone B-L43 857.090

Product Name Anti-Human IL1RAP Monoclonal Antibody, Azide Free Clone B-L43

Catalog Number 857.090

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, FC

Immunogen Recombinant human IL1RAP

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognizes natural and recombinant human IL1RAP

References **Trad, R. et al., J Immunother Cancer. 2022; 10(7): e004222.**

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells