

Product Datasheet

Anti-Human CD64 Monoclonal Antibody, Azide Free Clone B-T44 857.120

Product Name Anti-Human CD64 Monoclonal Antibody, Azide Free Clone B-T44

Catalog Number 857.120

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications FC

Immunogen Recombinant Human CD64

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. DO NOT FREEZE.

Isotype Mouse IgG2b Kappa light chain

Clonality Monoclonal

Specificity Recognize natural and recombinant human CD64

Hybridoma Myeloma X63/AG.8653 x Balb/c lymph node cells