

Product Datasheet

Anti-Human IL-8 Monoclonal Antibody, Azide Free Clone B-B49 857.170

Product Name	Anti-Human IL-8 Monoclonal Antibody, Azide Free Clone B-B49
Catalog Number	857.170
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human IL-8
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant antibody from a phage display library
Specificity	Recognises both natural and recombinant human IL-8