

Product Datasheet

Anti-Human FGF2 Monoclonal Antibody, Azide Free Clone B-G58 857.180

Product Name	Anti-Human FGF2 Monoclonal Antibody, Azide Free Clone B-G58
Catalog Number	857.180
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human FGF basic
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile filtered through 0.22µm. Carrier and preservative free.
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes natural and recombinant human FGF2
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells