

Product Datasheet

Anti-Human HLA Class II (DQ) Monoclonal Antibody, Azide Free Clone B-K27 857.210

| | |
|----------------------------|---|
| Product Name | Anti-Human HLA Class II (DQ) Monoclonal Antibody, Azide Free Clone B-K27 |
| Catalog Number | 857.210 |
| Description | Mouse monoclonal |
| Species/Host | Mouse |
| Conjugation | Azide free |
| Tested Applications | IP, FC |
| Immunogen | T4+ CTL clone HG-38 |
| Brand | Diaclone |
| Form/Appearance | Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free |
| Storage | Stable at +2-8°C for 12 months. For longer storage freeze aliquots. |
| Isotype | Mouse IgG2a Kappa light chain |
| Clonality | Monoclonal |
| Specificity | Recognises a 28 kDa and 33 kDa complex |
| Hybridoma | Myeloma X63/AG.8653 x Balb/c spleen cells |