

Product Datasheet

Anti-Human HLA Class I Monoclonal Antibody, Unconjugated Clone W6/32 857.223.020

Product Name	Anti-Human HLA Class I Monoclonal Antibody, Unconjugated Clone W6/32
Catalog Number	857.223.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IP, FC
Immunogen	Membrane of human tonsil cells
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises an antigenic determinant shared among products of the HLA-A, -B and -C loci
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells