

Product Datasheet

Anti-Human IL-28B Monoclonal Antibody, Azide Free Clone B-Z22 857.270

Product Name Anti-Human IL-28B Monoclonal Antibody, Azide Free Clone B-Z22

Catalog Number 857.270

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA

Immunogen Recombinant human IL-28B

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant human IL-28B

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells