

Product Datasheet

Anti-Human RANTES Monoclonal Antibody, Azide Free Clone B-C53 857.600

| | |
|----------------------------|--|
| Product Name | Anti-Human RANTES Monoclonal Antibody, Azide Free Clone B-C53 |
| Catalog Number | 857.600 |
| Description | Mouse monoclonal |
| Species/Host | Mouse |
| Conjugation | Azide Free |
| Tested Applications | ELISA |
| Immunogen | Recombinant human RANTES |
| Brand | Diaclone |
| Form/Appearance | Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free |
| Storage | Stable at +2-8°C for 12 months. For longer storage freeze aliquots. |
| Isotype | Mouse IgG2a Kappa light chain |
| Clonality | Monoclonal |
| Specificity | Recognizes both natural and recombinant human RANTES |
| Hybridoma | Myeloma X63/AG.8653 x Balb/c node cells |