

Product Datasheet

Anti-h Transferrin 3104 SPTN-5 100525

Product Name	Anti-h Transferrin 3104 SPTN-5
Catalog Number	100525
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	KA= 3.8×10^9 1/M; KD= 2.6×10^{-10} M (=0.26 nM)
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.7-7.4
Cross Reactivity	N/D
Specificity	Antibody recognizes human transferrin
Shelf Life	Unspecified