

Product Datasheet

Anti-Human CD146 Monoclonal Antibody, Azide Free Clone B-T46 857.630

Product Name Anti-Human CD146 Monoclonal Antibody, Azide Free Clone B-T46

Catalog Number 857.630

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications WB, FC

Immunogen Recombinant human CD146

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. DO NOT FREEZE.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognizes both recombinant and natural human CD146

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells