

Product Datasheet

Anti-Human CD146 Monoclonal Antibody, Unconjugated Clone B-T46 857.633.020

| | |
|----------------------------|---|
| Product Name | Anti-Human CD146 Monoclonal Antibody, Unconjugated Clone B-T46 |
| Catalog Number | 857.633.020 |
| Description | Mouse monoclonal |
| Species/Host | Mouse |
| Conjugation | Purified |
| Tested Applications | FC |
| Immunogen | Recombinant human CD146 |
| Brand | Diaclone |
| Form/Appearance | Phosphate-buffered saline with 1% BSA and 0.09% sodium azide |
| Storage | Stable at +2-8°C for 12 months. DO NOT FREEZE. |
| Note | Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood |
| Isotype | Mouse IgG1 Kappa light chain |
| Clonality | Monoclonal |
| Specificity | Recognizes both recombinant and natural human CD146 |
| Hybridoma | Myeloma X63/AG.8653 x Balb/c node cells |