

Product Datasheet

Anti-Human CD9 Monoclonal Antibody, Azide Free Clone TS9 857.750

Product Name	Anti-Human CD9 Monoclonal Antibody, Azide Free Clone TS9
Catalog Number	857.750
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	IP, WB, FC
Immunogen	HeLa cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the p24 antigen, a 24 kda protein
Hybridoma	Myeloma X63/AG. x Balb/c spleen cells