

Product Datasheet

Anti-Human CD53 Monoclonal Antibody, Azide Free Clone TS53 857.760

Product Name	Anti-Human CD53 Monoclonal Antibody, Azide Free Clone TS53
Catalog Number	857.760
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	IP, WB, FC
Immunogen	Jurkat cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the MEM-53 antigen, a 30-50 kda (smear) protein
Hybridoma	Myeloma X63/AG8. x Balb/c spleen cells