

Product Datasheet

Anti-Human CD82 Monoclonal Antibody, Azide Free Clone TS82b 857.790

Product Name Anti-Human CD82 Monoclonal Antibody, Azide Free Clone TS82b

Catalog Number 857.790

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications IP, WB, FC

Immunogen Jurkat and HEL cell lines

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Carrier and preservative free

Storage Stable at +2-8°C for 12 months. DO NOT FREEZE.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the CD82 antigen, a 60 kDa protein

Hybridoma Myeloma X63/AG8. x Balb/c spleen cells