

## Product Datasheet

### Anti-Human ADAM-10 Monoclonal Antibody, Azide Free Clone 11G2 857.800

<b>Product Name</b>	Anti-Human ADAM-10 Monoclonal Antibody, Azide Free Clone 11G2
<b>Catalog Number</b>	857.800
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Azide free
<b>Tested Applications</b>	IP, WB, FC
<b>Immunogen</b>	Jurkat cell line
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
<b>Storage</b>	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
<b>Isotype</b>	Mouse IgG1 Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Recognises the ADAM-10 protein, either the 70 kDa form or the 97 kDa precursor form
<b>Hybridoma</b>	Myeloma X63/AG8. x Balb/c spleen cells