

Product Datasheet

Anti-Human CD105 Monoclonal Antibody, Azide Free Clone B-A55 857.810

Product Name Anti-Human CD105 Monoclonal Antibody, Azide Free Clone B-A55

Catalog Number 857.810

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide Free

Tested Applications ELISA, WB, FC

Immunogen Recombinant human CD105

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant human CD105

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells