

## Product Datasheet

### Anti-Human CD105 Monoclonal Antibody, Azide Free Clone B-H45 857.820

**Product Name** Anti-Human CD105 Monoclonal Antibody, Azide Free Clone B-H45

**Catalog Number** 857.820

**Description** Mouse monoclonal

**Species/Host** Mouse

**Conjugation** Azide Free

**Tested Applications** ELISA

**Immunogen** Recombinant human CD105

**Brand** Diaclone

**Form/Appearance** Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

**Storage** Stable at +2-8°C for 12 months. DO NOT FREEZE.

**Isotype** Mouse IgG1 Kappa light chain

**Clonality** Monoclonal

**Specificity** Recognises both natural and recombinant human CD105

**Hybridoma** Myeloma X63/AG.8653 x Balb/c node cells