

Product Datasheet

Rat Anti-Human IL-5 Monoclonal Antibody, Azide Free Clone B-Z25 857.850

Product Name Rat Anti-Human IL-5 Monoclonal Antibody, Azide Free Clone B-Z25

Catalog Number 857.850

Description Mouse monoclonal

Species/Host Rat

Conjugation Azide free

Tested Applications ELISA

Immunogen COS-expressed recombinant hIL-5

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Rat IgG2a Kappa light chain

Clonality Monoclonal

Specificity Human IL-5