

Product Datasheet

Anti-Human TIM-3 Monoclonal Antibody, Azide Free Clone B-K48 857.860

Product Name	Anti-Human TIM-3 Monoclonal Antibody, Azide Free Clone B-K48
Catalog Number	857.860
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, FC
Immunogen	Recombinant Human TIM-3/Fc
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognise natural (both soluble and membrane forms) and recombinant human TIM-3
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells