

Product Datasheet

Anti-Human IFN alpha 2a Monoclonal Antibody, Azide Free Clone B-Z38 857.900

Product Name	Anti-Human IFN alpha 2a Monoclonal Antibody, Azide Free Clone B-Z38
Catalog Number	857.900
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human IFN alpha 2a
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes both natural and recombinant human interferon alpha 2a
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells