

Product Datasheet

Anti-Human IFN alpha 2a Monoclonal Antibody, Azide Free Clone B-N42 857.910

| | |
|----------------------------|--|
| Product Name | Anti-Human IFN alpha 2a Monoclonal Antibody, Azide Free Clone B-N42 |
| Catalog Number | 857.910 |
| Description | Mouse monoclonal |
| Species/Host | Mouse |
| Conjugation | Azide free |
| Tested Applications | ELISA |
| Immunogen | Recombinant human IFN alpha 2a |
| Brand | Diaclone |
| Form/Appearance | Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free. |
| Storage | Stable at +2-8°C for 12 months. DO NOT FREEZE. |
| Isotype | Mouse IgG1 Kappa light chain |
| Clonality | Monoclonal |
| Specificity | Recognizes both natural and recombinant human interferon alpha 2a |
| Hybridoma | Myeloma X63/AG.8653 x Balb/c node cells |