

Product Datasheet

Anti-Human Galectin-9 Monoclonal Antibody, Unconjugated Clone B-T33 857.943.020

Product Name	Anti-Human Galectin-9 Monoclonal Antibody, Unconjugated Clone B-T33
Catalog Number	857.943.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant human Galectin-9
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human Galectin-9 protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells