

Product Datasheet

Anti-Human ST2 Monoclonal Antibody, Unconjugated Clone B-L42 857.953.020

Product Name Anti-Human ST2 Monoclonal Antibody, Unconjugated Clone B-L42

Catalog Number 857.953.020

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested Applications FC

Immunogen Recombinant human ST2

Brand Diaclone

Form/Appearance Phosphate-buffered saline with 1% BSA and 0.09% sodium azide

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Note Use 10 µl to label 10⁶ cells or 100 µl of whole blood

Isotype Mouse IgG2b Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant human ST2

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells