

Product Datasheet

Anti-Human LL-37 Monoclonal Antibody, Azide Free Clone B-T37 857.960

Product Name Anti-Human LL-37 Monoclonal Antibody, Azide Free Clone B-T37

Catalog Number 857.960

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, FC

Immunogen Human antimicrobial peptide LL-37

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. DO NOT FREEZE.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant human LL-37

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells