

Product Datasheet

Anti-Human LL-37 Monoclonal Antibody, Azide Free Clone B-F57 857.980

Product Name	Anti-Human LL-37 Monoclonal Antibody, Azide Free Clone B-F57
Catalog Number	857.980
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, WB
Immunogen	Human antimicrobial peptide LL-37
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human LL-37
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells