

Product Datasheet

Accessory pack for ELISA set 858.000.005

Product Name	Accessory pack for ELISA set
Catalog Number	858.000.005
Description	<p>Material for optimal assay performance (quantity provided for 5x96 tests)</p> <p>Contents:</p> <p>5 x 96 wells strip plates and plate storage material, Coating Buffer 10X 1 vial 5ml, Wash Buffer 200X - 1 vial 25ml, Blocking Buffer - 5 vials 10ml, Standard Dilution Buffer 10X - 1 vial 12.5ml, Reconstitution Buffer 10X - 1 vial 1 ml, Secondary Antibody Dilution Buffer - 1 vial 25 ml, HRP Dilution Buffer - 2 vials 25 ml, Stop reagent - 5 vials 10 ml</p>
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C