

Product Datasheet

Diaclone Lysing Buffer 858.040.125

| | |
|------------------------|---|
| Product Name | Diaclone Lysing Buffer |
| Catalog Number | 858.040.125 |
| Description | Erythrocyte lysing solution for preparation of leukocytes in whole blood for their analysis by flow cytometry. |
| Brand | Diaclone |
| Form/Appearance | Lyophilised |
| Storage | Store at room temperature. After reconstitution store at +2°C to 8°C. |
| Note | Reconstitute with 250 mL of distilled water, check and adjust the pH to 7.3. After reconstitution, the lysing solution is stable for 12 months. |
| Specificity | Erythrocyte lysing solution |