

Product Datasheet

Murine GM-CSF ELISA Kit 860.060

Product Name	Murine GM-CSF ELISA Kit
Catalog Number	860.060
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	15.6-500 pg/ml
Sensitivity	9 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes both natural and recombinant murine GM-CSF
Sample Types	Serum, Cell culture supernatant