

Product Datasheet

Murine IL-17A ELISA Kit 860.070

Product Name	Murine IL-17A ELISA Kit
Catalog Number	860.070
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	15.6-500 pg/ml
Sensitivity	4.6pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine IL-17A
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable