

## Product Datasheet

### Murine CD138 ELISA Kit 860.090

<b>Product Name</b>	Murine CD138 ELISA Kit
<b>Catalog Number</b>	860.090
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	0.5-16 ng/ml
<b>Sensitivity</b>	0.35 ng/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	Syndecan-1
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other molecules
<b>Specificity</b>	Recognizes recombinant and natural murine soluble Syndecan-1
<b>Sample Types</b>	Cell culture supernatant, serum, plasma, buffered solution