

Product Datasheet

Murine CD138 ELISA Set 861.060

Product Name	Murine CD138 ELISA Set
Catalog Number	861.060
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	0.5-16 ng/ml
Sensitivity	0.35 ng/ml
Tested Applications	ELISA
Alternative Names	Syndecan-1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other murine soluble molecule
Specificity	Recognizes recombinant and natural murine soluble Syndecan-1
Sample Types	Cell culture supernatant, serum, plasma