

## Product Datasheet

### Human IL-10 ELISpot Pair 869.030.010

<b>Product Name</b>	Human IL-10 ELISpot Pair
<b>Catalog Number</b>	869.030.010
<b>Description</b>	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	NULL
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10.
<b>Specificity</b>	Recognizes natural human IL-10