

Product Datasheet

Human TNF alpha ELISpot Pair 869.050.010

| | |
|----------------------------|--|
| Product Name | Human TNF alpha ELISpot Pair |
| Catalog Number | 869.050.010 |
| Description | Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples. |
| Tested Applications | ELISpot |
| Alternative Names | TNFa |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with soluble TNFR I and TNFR II. |
| Specificity | Recognizes natural human TNF-a |