

Product Datasheet

Human IFN gamma ELISpot Pair 869.060.010

Product Name	Human IFN gamma ELISpot Pair
Catalog Number	869.060.010
Description	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g.
Specificity	Recognizes natural human IFN-g