

Product Datasheet

Human IL-1 beta ELISpot Pair 869.100.010

Product Name	Human IL-1 beta ELISpot Pair
Catalog Number	869.100.010
Description	Diaclone ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Tested Applications	ELISpot
Alternative Names	IL-1b, IL-1beta
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes natural (pro and mature) human IL-1b