

## Product Datasheet

### Rat IFN gamma ELISA Set 872.000

<b>Product Name</b>	Rat IFN gamma ELISA Set
<b>Catalog Number</b>	872.000
<b>Description</b>	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
<b>Range</b>	31.25-1000 pg/ml
<b>Sensitivity</b>	10 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	IFN-g, IFN-gamma
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other rat cytokines
<b>Specificity</b>	Recognizes both natural and recombinant rat IFN-g
<b>Sample Types</b>	Serum, Cell culture supernatant