

Product Datasheet

Human IL-29 ELISA Kit 873.050

Product Name	Human IL-29 ELISA Kit
Catalog Number	873.050
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	250 pg/ml - 7.8 pg/ml
Sensitivity	5.9 pg/ml
Tested Applications	ELISA
Alternative Names	Interferon lambda 1, ZCYTO21
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity observed with IFNa2, IL-28A, IL-28B, IFNg, IL-10, IL-24, IP-10, IL-31, IL-34 and GM-CSF
Specificity	Recognizes both natural and recombinant human IL-29
Sample Types	Cell culture supernatant, serum or plasma tested, other body fluids suitable