

Product Datasheet

Human IL-34 ELISA Kit 873.060

Product Name	Human IL-34 ELISA Kit
Catalog Number	873.060
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	5000 pg/ml - 156.25 pg/ml
Sensitivity	39 pg/ml
Tested Applications	ELISA
Alternative Names	Protein C16orf77
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity observed with CD31, IL-10, IP-10, IL-31, IL-1a, IL-1b, GM-CSF, M-CSF, IL-29, MCP-1 and murine IL-34
Specificity	Recognizes both natural and recombinant IL-34
Sample Types	Cell culture supernatant, serum or plasma tested, other body fluids suitable