

## Product Datasheet

### Human Angiopoietin-2 ELISA Kit 873.080

<b>Product Name</b>	Human Angiopoietin-2 ELISA Kit
<b>Catalog Number</b>	873.080
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	125-4000 pg/ml
<b>Sensitivity</b>	25 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	ANG-2,ANGPT-2
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with VEGF-A, GM-CSF, IL-34, IL-29 IFNg, IL-2, IL-10, CD138, MICB, MCP-1
<b>Specificity</b>	Recognizes natural and recombinant Angiopoietin-2
<b>Sample Types</b>	Serum, Cell culture supernatant