

Product Datasheet

Human MICB ELISA Kit 873.090

Product Name	Human MICB ELISA Kit
Catalog Number	873.090
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	12.3 pg/ml
Tested Applications	ELISA
Alternative Names	PERB11.2
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with MICA, Ang-2, IL-1b, IL-34, IL-10, IL-31, IFNg, GM-CSF, VEGF-165 and IL-2
Specificity	Recognizes both natural and recombinant human MICB
Sample Types	Serum/Plasma/Cell culture supernatants