

Product Datasheet

Human Galectin-9 ELISA Kit 873.100

Product Name	Human Galectin-9 ELISA Kit
Catalog Number	873.100
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	15.6-500 pg/ml
Sensitivity	8 pg/ml
Tested Applications	ELISA
Alternative Names	Ecalectin, HUAT, LGALS9
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with GM-CSF, MICB, IL-34, IFN γ , IL-29, Galectin-1, Galectin-3, Galectin-8, Galectin-10
Specificity	Recognizes both natural and recombinant human Galectin-9
Sample Types	Serum/Plasma/Cell culture supernatants