

Product Datasheet

Human IFN γ / IL-4 Dual Fluorospot Set 874.002.005

| | |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Human IFN γ / IL-4 Dual Fluorospot Set |
| Catalog Number | 874.002.005 |
| Description | Diaclone Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection antibody for cytokine 1, anti-FITC antibody green fluorescence conjugate, capture antibody for cytokine 2, biotinylated detection antibody for cytokine 2, streptavidin-phycoerythrin conjugate, BSA. Kit can be sold with or without plates. |
| Tested Applications | Fluorospot |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | Cross reactivity with simian IFN-g |
| Specificity | Recognizes natural human IL-4 and IFN-g |