
Product Datasheet

Human IFN γ / IL-2 Dual ELISpot Set 874.040.005

Product Name	Human IFN γ / IL-2 Dual ELISpot Set
Catalog Number	874.040.005
Description	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, and blocking reagent.
Tested Applications	ELISpot
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	Cross reactivity with simian IFN- γ
Specificity	Recognizes both natural and recombinant human IFN- γ and IL-2