

## Product Datasheet

### Human IFN $\gamma$ / IL-2 Dual Fluorospot Set 874.042.010

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Human IFN $\gamma$ / IL-2 Dual Fluorospot Set  |
| <b>Catalog Number</b>      | 874.042.010  |
| <b>Description</b>         | Diaclone Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection antibody for cytokine 1, anti-FITC antibody green fluorescence conjugate, capture antibody for cytokine 2, biotinylated detection antibody for cytokine 2, streptavidin-phycoerythrin conjugate, BSA. Kit can be sold with or without plates. |
| <b>Tested Applications</b> | Fluorospot   |
| <b>Brand</b>               | Diaclone   |
| <b>Storage</b>             | Store at +4°C  |
| <b>Cross Reactivity</b>    | Cross reactivity with simian IFN-g   |
| <b>Specificity</b>         | Recognizes both natural and recombinant human IFN-g and IL-2   |