

Product Datasheet

Anti-Human IL-10 Detection Antibody, Biotin Conjugated Clone B-T10 879.060.002

Product Name	Anti-Human IL-10 Detection Antibody, Biotin Conjugated Clone B-T10
Catalog Number	879.060.002
Description	This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA, ELISpot
Immunogen	Recombinant human IL-10
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide
Storage	Stable at +2-8°C for 12 months.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA/ELISpot: 0,05-0,5µg/ml
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-10
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells