

Product Datasheet

Anti-Human IL-12 p35 Capture Antibody (IL-12 p70 ELISA), Azide Free Clone B-T21 879.070.001

Product Name	Anti-Human IL-12 p35 Capture Antibody (IL-12 p70 ELISA), Azide Free Clone B-T21
Catalog Number	879.070.001
Description	This Antibody can be used as Capture in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, ELISpot
Immunogen	Recombinant human IL-12
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,3-1µg/ml ELISpot: 5-10µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-12 p35 + p70
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells