

## Product Datasheet

### Anti-Human IL-13 Detection Antibody, Biotin Conjugated Clone B-P6 879.080.002

<b>Product Name</b>	Anti-Human IL-13 Detection Antibody, Biotin Conjugated Clone B-P6
<b>Catalog Number</b>	879.080.002
<b>Description</b>	This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay.
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, ELISpot
<b>Immunogen</b>	Recombinant human IL-13
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide
<b>Storage</b>	Stable at +2-8°C for 12 months.
<b>Note</b>	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA/ELISpot: 0,05-0,5µg/ml
<b>Isotype</b>	Mouse IgG1 Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Recognises both natural and recombinant human IL-13
<b>Hybridoma</b>	Myeloma X63/AG.8653 x Balb/c spleen cells