

Product Datasheet

Anti-Human TNF α Detection Antibody, Biotin Conjugated Clone B-C7 879.090.002

Product Name	Anti-Human TNF α Detection Antibody, Biotin Conjugated Clone B-C7
Catalog Number	879.090.002
Description	This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA, ELISpot
Immunogen	Recombinant human TNF-a
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide
Storage	Stable at +2-8°C for 12 months.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA/ELISpot: 0,05-0,5 μ g/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human TNF-a
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells