

Product Datasheet

Anti-Human IL-10 Capture Antibody, Azide Free Clone B-S10 879.100.001

Product Name Anti-Human IL-10 Capture Antibody, Azide Free Clone B-S10

Catalog Number 879.100.001

Description This Antibody can be used as Capture in a human sandwich Immunoassay.

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, ELISpot

Immunogen Recombinant human IL-10

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Note The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 1-5µg/ml ELISpot: 5-10µg/ml

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant human IL-10

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells