

## Product Datasheet

### Anti-Human IL-34 Capture Antibody, Azide Free Clone B-Z34 879.120.001

**Product Name** Anti-Human IL-34 Capture Antibody, Azide Free Clone B-Z34

**Catalog Number** 879.120.001

**Description** This Antibody can be used as Capture in a human sandwich Immunoassay.

**Species/Host** Mouse

**Conjugation** Azide free

**Tested Applications** ELISA

**Immunogen** Recombinant human IL-34

**Brand** Diaclone

**Form/Appearance** Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

**Storage** Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

**Note** The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 1-10µg/ml

**Isotype** Mouse IgG1 Kappa light chain

**Clonality** Monoclonal

**Specificity** Recognises both natural and recombinant human IL-34

**Hybridoma** Myeloma X63/AG.8653 x Balb/c node cells